

ACUSERA 

# Blood Gas Controls

---

Blood Gas Product Range		
Product Description	Size	Cat. No.
Blood Gas Control Level 1	30 x 1.8 ml	BG5001
Blood Gas Control Level 2	30 x 1.8 ml	BG5002
Blood Gas Control Level 3	30 x 1.8 ml	BG5003

## Blood Gas Control

Analytes			
pCO <sub>2</sub> pH pO <sub>2</sub>	Total CO <sub>2</sub> Calcium Chloride	Glucose Lactate Potassium	Sodium

Liquid ready-to-use blood gas control designed for use in the routine monitoring of the 10 parameters listed above. Method specific target values are provided for the most common instruments.

- **Aqueous material**
- **Stable to expiry date at 2°C to 8°C**
- **Once opened controls should be analysed within 1 minute for pH and blood gas parameters; for electrolyte measurements, the control should be analysed within 1 hour of opening**

Description	Size	Cat. No.
Blood Gas Control Level 1	30 x 1.8 ml	BG5001
Blood Gas Control Level 2	30 x 1.8 ml	BG5002
Blood Gas Control Level 3	30 x 1.8 ml	BG5003