

RANDOX

COVID-19 RAPID TEST

Results in 39 Minutes

Our SARS-CoV-2 (COVID-19) rapid test is a CE marked real time PCR cartridge, providing clear and concise results in a timely manner, direct at the point of care. This enables individuals to take the recommended safety precautions without delay. The SARS-CoV-2 rapid test targets the E gene sequence for COVID-19 detection in just 39 minutes with minimal workflow required, suitable for laboratory and non-laboratory settings.

Detection Method: Real-Time PCR


Result Time: 39 Minutes

Sample Volume: 300µl Clinical Sample

Sample Type: Nasopharyngeal or Oropharyngeal Swab

VIRUS

SARS-CoV-2 (COVID-19)



RANDOX

VIVALYTIC

Cartridge Based Testing Platform


The new SARS-CoV-2 rapid test will be conducted on Vivalytic, a point of care platform brought to the market by Randox Laboratories and Bosch. The Vivalytic system is a fully-automated, cartridge-based platform capable of both Hi-Plex and Lo-Plex infectious disease testing. Each easy to use cartridge contains all necessary reagents, is fully-sealed to minimise risk and can be conveniently stored at room temperature.

Vivalytic Key Features:

- Cartridge based platform
- Fully-automated
- Suitable for both non-laboratory & laboratory settings
- Small footprint device
- Rapid turnaround from sample entry to results
- Available as stand-alone or modular with up to 8 platforms with Vivalytic Up

VIVALYTIC WORKFLOW

4 Easy Steps for Optimised Workflow



To find out more email:
marketing@randox.com